

## What is ClearV™?

ClearV™ with its unique integrated vein illumination technology allows your physician to precisely treat vessels, broken capillaries, spider veins, vascular lesions, and discoloration of your skin that you can and can't see!

## What will ClearV™ do?

Your physician will use ClearV to restore your skin tone and color. Broken capillaries and other visible vascular structures will be gently eliminated from the upper layers of your skin surface. Unwanted vessels will be heated by the laser energy and become less noticeable.

## What conditions can be treated?

- Broken Capillaries
- Spider Veins
- Telangiectasia
- Facial Veins
- Rosacea



ClearV™

LOOKING TO  
SHOW OFF  
YOUR LEGS

WITHOUT  
UNSIGHTLY  
VEINS

 SCITON®  
BECAUSE RESULTS MATTER

925 Commercial Street  
Palo Alto, California 94303  
(888) 646-6999 • [info@sciton.com](mailto:info@sciton.com)

©2019 Sciton, Inc. All rights reserved. Sciton is a registered trademark and ClearV and SVI are trademarks of Sciton, Inc.  
Printed in USA

2600-013-01 Rev A

# ClearV™ Results



ClearV™ | 4 Month Post 1 tx | courtesy of Farzad Tabibzadeh, MD



ClearV™ | 1 Month Post 1 tx | courtesy of Brian V. Heil, MD



ClearV™ | Immediately After 1 tx | courtesy of Dana Miriciouis, MD



ClearV™ | 3 Months Post 2 tx | courtesy of Arkansas Laser Solutions

## TREAT THE VEINS YOU CAN AND CAN'T SEE!

**ClearV™ Subsurface Vein Illumination (SVI™)** – Specially designed to enhance the visibility to see subsurface and feeder veins to treat more thoroughly and effectively for enhanced treatment precision and patient comfort.



## FREQUENTLY ASKED QUESTIONS

### What areas can be treated?

All areas can be treated, although the more popular areas are the legs and face. Discuss your needs with your physician.

### How does vein illumination work?

The ClearV handpiece is uniquely designed to visualize sub-surface veins with yellow or orange light which are integrated into the handpiece. This provides the treating physician with enhanced visibility, allowing for a more thorough and effective treatment.

### What will happen during the procedure?

Your eyes will be protected with safety shields or glasses. Laser energy will be delivered precisely to the area of your skin being treated. You may feel a brief warm or hot sensation as the laser pulses are absorbed by the unwanted vessels. Your procedure may take a few minutes to half an hour depending on the size of the area being treated.

### Will it be uncomfortable?

There should be little to no discomfort, however it varies from patient to patient and depends on the treatment. There is no need for anesthesia and your physician will discuss your options to properly manage any discomfort from the treatment. ClearV has integrated cooling which will improve your experience. There should be no prolonged post treatment discomfort.

### What should I expect after treatment?

Patient responses can vary. Redness and swelling around the treated vessels are usually noted within a few minutes and typically have a “cat scratch” appearance. A sun burn sensation and a slight feeling of soreness are also normal and expected. These reactions tend to subside within 24-48 hours after the treatment.

After 48 hours, the skin over the treated veins will look slightly red, purple, or bruised. Sometimes the veins will be more visible than they were prior to treatment. This is normal. Within the first two to three weeks, the redness and purple color will begin to fade. Each week, the color will become lighter and the veins will become less noticeable.

### What aftercare will I need?

It is important to follow your physician's instructions. You may be given some skin care products and instructions on how to use them. Strictly avoid any sun exposure to the treated area for a minimum of 7 days after the procedure. The treated area is more prone to sunburn and pigmentation change. Keep the area covered and use a sun block with a protection factor of 30+ and reapply the sun block every two hours.

### How long will it take to recover and how many treatments do I need?

Your healing time will depend upon your actual treatment. In most cases, there is virtually no downtime and you may resume your daily activities immediately. Subsequent treatments are based upon your clinician's recommendation and are typically 4-6 weeks apart. The number of treatments required will depend upon the number of vessels present, the character of the vessels, and your body's ability to heal.